

BLOOD COLLECTING SYSTEMS

Reference	Description	Size (mm)	Additive	Cap color/ color	Shelf life (months)	Pack/ Carton
	VACUSERA BLOOD COLLECTION TUBES					
234104	VACUSERA 4 ml No Additive	13x75	No additive	WHITE	18 months	100/1200
235106	VACUSERA 6 ml No Additive	13x100	No additive	WHITE	18 months	100/1200
236109	VACUSERA 9 ml No Additive	16x100	No additive	WHITE	18 months	100/1200
234202	VACUSERA 2 ml Serum Clot Activator	13x75	Clot Activator	RED	18 months	100/1200
234203	VACUSERA 3 ml Serum Clot Activator	13x75	Clot Activator	RED	18 months	100/1200
234204	VACUSERA 4 ml Serum Clot Activator	13x75	Clot Activator	RED	18 months	100/1200
235206	VACUSERA 6 ml Serum Clot Activator	13x100	Clot Activator	RED	18 months	100/1200
236209	VACUSERA 9 ml Serum Clot Activator	16x100	Clot Activator	RED	18 months	100/1200
236200	VACUSERA 10 ml Serum Clot Activator	16x100	Clot Activator	RED	18 months	100/1200
234303	VACUSERA 3.5 ml Serum Gel and Clot Activator	13x75	GEL + Clot Activator	ORANGE	18 months	100/1200
235305	VACUSERA 5 ml Serum Gel and Clot Activator	13x100	GEL + Clot Activator	ORANGE	18 months	100/1200
235315	VACUSERA 5 ml Serum Gel and Thrombin	13x100	GEL + Thrombin	ORANGE	18 months	100/1200
236308	VACUSERA 8 ml Serum Gel and Clot Activator	16x100	GEL + Clot Activator	ORANGE	18 months	100/1200
234502	VACUSERA 1.8 ml Coagulation %3.2 Sodium Citrate - Double Tube	13x75	9NC Sodium Citrate 3.2%(0.109M)	BLUE	9 months	100/1200
234503	VACUSERA 2.7 ml Coagulation %3.2 Sodium Citrate - Double Tube	13x75	9NC Sodium Citrate 3.2%(0.109M)	BLUE	9 months	100/1200
234504	VACUSERA 3,6 ml Coagulation %3.2 Sodium Citrate	13x75	9NC Sodium Citrate 3.2%(0.109M)	BLUE	4 months	100/1200
234512	VACUSERA 1.8 ml Coagulation %3.8 Sodium Citrate - Double Tube	13x75	9NC Sodium Citrate 3.8%(0.129M)	BLUE	9 months	100/1200
234513	VACUSERA 2.7 ml Coagulation %3.8 Sodium Citrate - Double Tube	13x75	9NC Sodium Citrate 3.8%(0.129M)	BLUE	9 months	100/1200
234514	VACUSERA 3,6 ml Coagulation %3.8 Sodium Citrate	13x75	9NC Sodium Citrate 3.8%(0.129M)	BLUE	4 months	100/1200
234523	VACUSERA 2,7 ml CTAD (SodiumCitrate %3,2)	13x75	Buffered Sodium Citrate	BLUE	9 months	100/1200
235006	VACUSERA 6 ml Trace Element(Clots Activator)	13x100	Clot Activator	NAVY BLUE	18 months	100/1200
236726	VACUSERA 6ml ACD-A blood grouping tubes	13x100	ACD-A	YELLOW	12 months	100/1200
236729	VACUSERA 9ml ACD-A blood grouping tubes	16x100	ACD-A	YELLOW	12 months	100/1200
236736	VACUSERA 6ml ACD-B blood grouping tubes	13x100	ACD-B	YELLOW	12 months	100/1200

BLOOD COLLECTING SYSTEMS

236739	VACUSERA 9ml ACD-B blood grouping tubes	16x100	ACD-B	YELLOW	12 months	100/1200
234602	VACUSERA 2 ml K2EDTA	13x75	K2E K2EDTA	VIOLET	18 months	100/1200
234603	VACUSERA 3 ml K2EDTA	13x75	K2E K2EDTA	VIOLET	18 months	100/1200
234604	VACUSERA 4 ml K2EDTA	13x75	K2E K2EDTA	VIOLET	18 months	100/1200
234605	VACUSERA 4.5 ml K2EDTA	13x75	K2E K2EDTA	VIOLET	18 months	100/1200
235606	VACUSERA 6 ml K2EDTA	13x100	K2E K2EDTA	VIOLET	18 months	100/1200
236609	VACUSERA 9 ml K2EDTA	16x100	K2E K2EDTA	VIOLET	18 months	100/1200
234623	VACUSERA 3 ml K2EDTA+ GEL	13x75	Gel + K2E K2EDTA	VIOLET	18 months	100/1200
235625	VACUSERA 5 ml K2EDTA+ GEL	13x100	Gel + K2E K2EDTA	VIOLET	18 months	100/1200
236628	VACUSERA 8 ml K2EDTA+ GEL	16x100	Gel + K2E K2EDTA	VIOLET	18 months	100/1200
234612	VACUSERA 2 ml K3EDTA	13x75	K3E K3EDTA	VIOLET	18 months	100/1200
234613	VACUSERA 3 ml K3EDTA	13x75	K3E K3EDTA	VIOLET	18 months	100/1200
234614	VACUSERA 4 ml K3EDTA	13x75	K3E K3EDTA	VIOLET	18 months	100/1200
234615	VACUSERA 4.5 ml K3EDTA	13x75	K3E K3EDTA	VIOLET	18 months	100/1200
235616	VACUSERA 6 ml K3EDTA	13x100	K3E K3EDTA	VIOLET	18 months	100/1200
236619	VACUSERA 9 ml K3EDTA	16x100	K3E K3EDTA	VIOLET	18 months	100/1200
235806	VACUSERA 6 ml Urine No Additive	13x100	NO ADDITIVE	YELLOW	18 months	100/1200
236800	VACUSERA 10 ml Urine No Additive	16x100	NO ADDITIVE	YELLOW	18 months	100/1200
235816	VACUSERA 6 ml Urine - Boric Acid	13x100	Boric Acid	YELLOW	18 months	100/1200
236810	VACUSERA 10 ml Urine - Boric Acid	16x100	Boric Acid	YELLOW	18 months	100/1200
235820	VACUSERA 10 ml Urine - Chlorhexidine	16x100	Chlorhexidine	YELLOW	18 months	100/1200
234402	VACUSERA 2 ml Li-Heparin	13x75	LH/Lithium Heparin	GREEN	18 months	100/1200
234404	VACUSERA 4 ml Li-Heparin	13x75	LH/Lithium Heparin	GREEN	18 months	100/1200
235406	VACUSERA 6 ml Li-Heparin	13x100	LH/Lithium Heparin	GREEN	18 months	100/1200
236409	VACUSERA 9 ml Li-Heparin	16x100	LH/Lithium Heparin	GREEN	18 months	100/1200
234424	VACUSERA 4 ml NH-Heparin	13x75	NH/Sodium Heparin	GREEN	18 months	100/1200
235426	VACUSERA 6 ml NH-Heparin	13x100	NH/Sodium Heparin	GREEN	18 months	100/1200
236429	VACUSERA 9 ml NH-Heparin	16x100	NH/Sodium Heparin	GREEN	18 months	100/1200
234413	VACUSERA 3 ml Li-Heparin with GEL	13x75	GEL + LH/Lithium Heparin	GREEN	18 months	100/1200
235415	VACUSERA 5 ml Li-Heparin with GEL	13x100	GEL + LH/Lithium Heparin	GREEN	18 months	100/1200

BLOOD COLLECTING SYSTEMS

236418	VACUSERA 8 ml Li-Heparin with GEL	16x100	GEL + LH/Lithium Heparin	GREEN	18 months	100/1200
234702	VACUSERA 2 ml Sodium Fluoride / K3EDTA	13x75	FE/Sodium Fluoride +EDTA K3	GRAY	18 months	100/1200
234704	VACUSERA 4 ml Sodium Fluoride / K3EDTA	13X75	FE/Sodium Fluoride +EDTA K3	GRAY	18 months	100/1200
234712	VACUSERA 2 ml Sodium Fluoride / Potassium Oxalate	13x75	FE/Sodium Fluoride +Potassium Oxalate	GRAY	18 months	100/1200
234714	VACUSERA 4 ml Sodium Fluoride / Potassium Oxalate	13X75	FE/Sodium Fluoride +Potassium Oxalate	GRAY	18 months	100/1200
234722	VACUSERA 2 ml Sodium Fluoride / Na2EDTA	13x75	FE/Sodium Fluoride + Na2EDTA	GRAY	18 months	100/1200
234724	VACUSERA 4 ml Sodium Fluoride / Na2EDTA	13X75	FE/Sodium Fluoride + Na2EDTA	GRAY	18 months	100/1200
	VACUSERA ESR- SEDIMENTATION SYSTEMS					
234902	VACUSERA 2 ml ESR Tube	13x75	4NC Sodium Citrate 3.2%	BLACK	12 months	100/1200
234900	VACUSERA Disposable ESR Pipette For 2 ml tubes				60 months	200/2000
234901	VACUSERA ESR Rack for 2 ml tube (open system)					1
233902	VACUSERA 1,6 ml ESR Tube - Glass	9*120	4NC Sodium Citrate 3.2%	BLACK	24 months	100/1000
233903	VACUSERA 1,28 ml ESR Tube - Glass	8*120	4NC Sodium Citrate 3.2%	BLACK	24 months	100/1000
233901	VACUSERA ESR Rack for 1,6 ml and 1,28 ml tube (closed system)					1
	VACUSERA BLOOD COLLECTION NEEDLES					
250018K	VACUSERA Blood Collection Needle 18 G x 1 1/2" (0,9x38mm)	18 G x1 1/2"	(1.2 x38mm)	PINK	60 months	100/5000
250020K	VACUSERA Blood Collection Needle 20 G x 1 1/2" (0,9x38mm)	20 G x1 1/2"	(0.9 x38mm)	YELLOW	60 months	100/5000
250021K	VACUSERA Blood Collection Needle 21 G x 1 1/2" (0,8x38mm)	21 G x1 1/2"	(0.8*38mm)	GREEN	60 months	100/5000
250022K	VACUSERA Blood Collection Needle 22 G x 1 1/2" (0,7x38mm)	22 G x1 1/2"	(0.7 x38mm)	BLACK	60 months	100/5000
255021K	VACUSERA Blood Collection Needle 21 G VISION	21 G x1 1/2"	(0.8*38mm)	GREEN	60 months	100/5000
255022K	VACUSERA Blood Collection Needle 22 G VISION	22 G x1 1/2"	(0.7 x38mm)	BLACK	60 months	100/5000
252020	VACUSERA Luer Adapter 20G	20 G			60 months	100/5000
251021K	VACUSERA Integrated Safety Needle 21 G	21G	21 G x 1 1/4" (0.8x32mm)	GREEN	60 months	48/2400
251022K	VACUSERA Integrated Safety Needle 22 G	22G	22 G x 1 1/4" (0.7x32mm)	BLACK	60 months	48/2400
260021K	Blood Collection Needle 21G (0,8x38mm) + Safety Holder	21G	21 G x 1 1/4" (0.8x32mm)	GREEN	60 months	25pcs/pack. , 500 case
260001	VACUSERA Standard Holder					200 /1000
260002	VACUSERA Safety Holder					50/800

BLOOD COLLECTING SYSTEMS

260005	VACUSERA Holder With Male Luer Lock					200/800
259001	VACUSERA Automatic Tourniquet				500	
259001-1	VACUSERA Automatic Tourniquet For Children - Green				500	
259001-2	VACUSERA Automatic Tourniquet For Children - Pink				500	
259050	VACUSERA Disposable Tourniquet				50	
	VACUSERA MINI CAPILLARY BLOOD COLLECTION SYSTEM					
238205	VACUSERA MINI 0.5 ml Serum Clot Activator		Clot Activator	RED	18 months	100/1000
238305	VACUSERA MINI 0.5 ml Gel and Clot Activator		GEL + Clot Activator	GOLD	18 months	100/1000
238605	VACUSERA MINI 0.5 ml K2EDTA		K2EDTA	VIOLET	18 months	100/1000
238615	VACUSERA MINI 0.5 ml K3EDTA		K3EDTA	VIOLET	18 months	100/1000
	VACUSERA SAFETY LANCET					
2393018	VACUSERA Safety Lancet 30G/1.8mm	30G/1.8mm		TURQUISE	36 months	100/1000 /5000
2392118	VACUSERA Safety Lancet 21G/1.8mm	21G/1.8mm		LIGHT BLUE	36 months	100/1000 /5000
2391818-BN	VACUSERA Safety Lancet - BLADE 18G/1.8mm	18G/1.8mm		PURPLE	36 months	100/1000 /5000
	VACUSERA SAFETY BLOOD COLLECTION SETS					
24121190K	VACUSERA Safety Blood Collection Set 21 G 3/4"x7"	0.8*19mm*180mm		GREEN	60 months	100/2000
24122190K	VACUSERA Safety Blood Collection Set 22 G 3/4"x7"	0.7*19mm*180mm		BLACK	60 months	100/2000
24123190K	VACUSERA Safety Blood Collection Set 23 G 3/4"x7"	0.6*19mm*180mm		BLUE	60 months	100/2000
24121300K	VACUSERA Safety Blood Collection Set with safety shield 21 G 3/4"x7"	0.8*19mm*300mm		GREEN	60 months	100/2000
24221190K	VACUSERA Safety Blood Collection Set with safety shield 21 G 3/4"x7"	0.8*19mm*180mm		GREEN	60 months	100/2000
24125190K	VACUSERA Safety Blood Collection Set 25 G 3/4"x7"	0.5*19mm*180mm		ORANGE	60 months	100/2000
	VACUSERA INTEGRATED URINE BEAKERS					
103023	VACUSERA URINE BEAKER WITH INTEGRATED HOLDER - NONSTERILE - WITHOUT BAG	100 ML		YELLOW	36 months	300
103026	VACUSERA URINE BEAKER WITH INTEGRATED HOLDER - STERILE - IN A SINGLE BAG	100 ML		YELLOW	36 months	300
103002	Urine container, red cap, non-sterile			RED	36 months	300
103003	Urine container, white cap, non-sterile			WHITE	36 months	300
236800	VACUSERA 10ml Urine Tube	16x100		YELLOW	18 months	100/2000
236820	VACUSERA 10ml Urine Tube- Chohexidine	16x100		YELLOW	18 months	100/2000
237101	VACUSERA PRP TUBE SET (BLUE, without gel, without accesories)	9 months				
237102	VACUSERA PRP TUBE SET WITH GEL (ORANGE, without accesories)	9 months				

BLOOD COLLECTING SYSTEMS

237111	VACUSERA PRP KIT SET (BLUE without gel, with Accesories)	9 months
237112A	VACUSERA PRP KIT SET WITH GEL (ORANGE, with accesories)	9 months
237202	VACUSERA PRF Tube	9 months
237201	VACUSERA PRF Tube	9 months



Contact

LABOVERITAS SIA

Reg. No 40203185225

VAT No LV40203185225

Address: Kleistu street 24, Riga LV-1067, LATVIA

<https://laboveritas.lv>

e-mail: info.laboveritas@gmail.com; andris@laboveritas.lv

mob.: +371 29161606

ISO 13485:2016